

**NEW!**  
Innovative  
one-piece top



The same  
tried and true  
frame as our  
576 model

## Roma® Overbed Tables

The new innovation behind our existing overbed table is the top. It has a one-piece, injection moulded polypropylene top without any voids for superior infection control (note that this is *not* a two-piece bonded product). It has smooth edges and internal corners for easy cleaning and reducing germ traps. This top has been developed and manufactured in New Zealand by Cubro in answer to the market requirements for a table that is easy clean, heat resistant and will have a greater product life in the demanding hospital environment.

**Ordering Code: 576M**

### INNOVATIVE FEATURES

#### One-piece top

One-piece, injection moulded, durable, polypropylene table top with CESA - antimicrobial.

#### CESA

An antimicrobial agent authorised by the FDA to be suitable for food contact application and offers:

- a complete spectrum against microorganism bacteria
- chemical and physical stability

#### Polypropylene

- does not absorb water like other plastics
- high heat resistance
- non-toxic (does not contain BPA found in polycarbonate plastics)
- stain resistant
- excellent resistance to cleaning chemicals
- durable material – won't delaminate (blister, chip or peel)

### STANDARD FEATURES

Powder coated base frame

Robust construction

Heat resistant to 160°

Smooth edges and internal corners for cleaning - no germ traps

Easy clean

Proven design with spring assisted, one hand adjustment

U-shape base fits under hospital beds

### TECHNICAL SPECIFICATIONS

|                   |                    |
|-------------------|--------------------|
| Top dimensions    | 828mm x 467mm      |
| Height            | 670mm - 1070mm     |
| Base height       | 100mm              |
| Inside base width | 740mm              |
| Castors           | 50mm (two braking) |