

# Compact oxygen concentrator



Easy to transport with smooth top and side handles



Increased patient safety with visible and audible alarms



## DESIGNED WITH ENHANCED QUALITY, DURABILITY AND SIMPLICITY IN MIND

The DeVilbiss 5 litre oxygen concentrator was designed with enhanced quality, durability and simplicity in mind.

The patented DeVilbiss TurnDown Technology minimizes wear on internal components, reduces power consumption by 15% or more and extends the life of the concentrator.

Model

Compact oxygen concentrator

Code

815

## FEATURES

- Oxygen Sensing Device is standard, ensuring patient safety and reliability for longer service intervals
- Audible and visual alarms indicate high/low pressure, low flow, low oxygen, power fail, and service required
- This helps to keep patients as safe as possible.
- Readily accessible patient controls for easier use
- Protected cannula fitting provides addition safety
- Recessed humidifier nook to prevent damage
- Pressure-compensated flow meter to ensure accurate flow display
- Pressure-relief valve and thermal protection on the compressor
- Convenient top and side handles for easy transport
- Flame-retardant cabinet

## SPECIFICATIONS

<p>W 342mm D 304mm H 622mm</p> <p>Overall dimensions</p>	<p>16.3kg</p> <p>Weight</p>
<p>0.5 to 5 LPM</p> <p>Delivery rate</p>	<p>The United States of America</p> <p>Assembled in</p>